

*DONJOY*®

GenuForce®

ELASTIC KNITTED KNEE SUPPORT  
WITH KNEECAP PAD



enovis™

# GenuForce®

## Elastic knitted knee support with kneecap pad

GenuForce's® contoured knee pad helps protect and support the patella, while its flexible bilateral stays help provide relief and stabilization for the knee joint.

### INTENDED USE/INDICATIONS:

- Support and compression following acute and chronic knee injury and mild knee osteoarthritis
- Proprioceptive and neuromuscular stimulation through compression
- Patellofemoral pain (knee pain)
- Mild knee instability
- Prevention of injury



### ANATOMICALLY CONTOURED SILICONE KNEE PAD

- for additional stability of the patella, massaging effect, and to improve proprioception




### FLEXIBLE BILATERAL STAYS

- help provide medial and lateral support



### ORDERING INFORMATION

DESCRIPTION	PART NO.	Ⓕ THIGH CIRCUMFERENCE	Ⓖ CALF CIRCUMFERENCE	SIZE
	82-1750-1	37 - 40 cm 14½ - 15¾"	25 - 28 cm 9¾ - 11"	1 = XS
	82-1750-2	40 - 43 cm 15¾ - 17"	28 - 31 cm 11 - 11¾"	2 = S
	82-1750-3	43 - 46 cm 17 - 18"	31 - 34 cm 11¾ - 13½"	3 = M
<b>GenuForce®</b>	82-1750-4	46 - 49 cm 18 - 19¼"	34 - 37 cm 13½ - 14½"	4 = L
	82-1750-5	49 - 52 cm 19¼ - 20½"	37 - 40 cm 14½ - 15¾"	5 = XL
	82-1750-6	52 - 55 cm 20½ - 21¾"	40 - 43 cm 15¾ - 17"	6 = XXL
	82-1750-7	55 - 58 cm 21¾ - 22¾"	43 - 46 cm 17 - 18"	7 = XXXL
	82-1751-5	52 - 55 cm 20½ - 21¾"	37 - 40 cm 14½ - 15¾"	5+
<b>GenuForce® Plus</b>	82-1751-6	55 - 58 cm 21¾ - 22¾"	40 - 43 cm 15¾ - 17"	6+
	82-1751-7	58 - 62 cm 22¾ - 24½"	43 - 46 cm 17 - 18"	7+
	82-1751-8	62 - 65 cm 24½ - 25½"	46 - 49 cm 18 - 19¼"	8+



### Ⓕ THIGH CIRCUMFERENCE

Measurement taken 15 cm above mid-patella

### Ⓖ CALF CIRCUMFERENCE

Measurement taken 5 cm below mid-patella

**enovis™**

Creating Better Together™

[enovisrehab.com.au](http://enovisrehab.com.au) / [enovisrehab.co.nz](http://enovisrehab.co.nz)

1/25 Frenchs Forest Road East, Frenchs Forest, NSW 2086

Individual results may vary. Neither Enovis, DJO, LLC or any of their subsidiaries dispense medical advice.

The contents of this document do not constitute medical, legal, or any other type of professional advice.

Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.

Copyright © 2023 DJO, LLC, a subsidiary of Enovis Corporation